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All raw materials used to produce our articles are conform with the most severe applicable Standards. Our Devices are:

A) suitable for foodstuff, according to EEC and International in vigour Directives; B) natural rubber, latex and CFC free; C) conform to EEC Directives and International Standards concerning heavy metals and phthalates limitations; D) biologically inert ( for example PE, PP and PC resulted non-cytotoxic). For any further information relevant to applicable Directives, Standards and Declarations of Conformity, you are kindly requested to contact our Quality Department (qualita@aptaca.com).

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### ACRYLONITRILE BUTADIENE STYRENE (ABS) ACRILONITRILE BUTADIENE STIRENE (ABS)

Compared to crystal polystyrene, that belongs to the same family, it's a more technical material and has a chemical resistance and temperature stability undoubtedly better. Products made of ABS are light and have a good rigidity grade. The most important mechanical properties are the high impact resistance and toughness also at low temperatures, hardness and scratch resistance. ABS is not autoclavable.

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### POLYSTYRENE (PS) POLISTIROLO O POLISTIRENE (PS)

It's a particularly transparent and bright material. It's rigid and hard but also rather fragile and subject to breakages. It is insensitive to moisture, steady at saline solutions, alkali, non oxidizing acids, aqueous chemical agents, but it has a very limited resistance to solvents. It's necessary to take into consideration that PS products, at ambient temperature, are particularly rigid and therefore they can crack or break in case they fall down from a laboratory bench or from higher heights. PS is not autoclavable.

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### POLYMETHYL METHACRYLATE (PMMA) POLIMETILMETACRILATO (PMMA)

PMMA is hard, rigid, quite fragile but extremely abrasion resistant (scratch resistant), it has high flexibility, it is stable to UV rays, but easily inflammable. PMMA has a good impact strength, higher than that of glass, but significantly lower than that of polycarbonate. Thanks to its optical properties, this plastic material is the most similar to glass. PMMA is not autoclavable.

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### POLYMETHYLPENTENE (PMP-TPX®) POLIMETILPENTENE (PMP - TPX®)

It's a very light resin, transparent (90% of transmittance of visual light) and very hard, with good mechanical properties. PMP resists to concentrated sulfuric acid, acetone, ethanol, mineral acids, strong and weak acids, strong and weak alkaline solutions, inorganic salts, aldehydes, alcohols, detergents, oils, fats and boiling water. It's, with a limited extent, steady to ketones. It has an optimal resistance to high temperatures (it can be repeatedly autoclaved both at +121 °C and at +150 °C). Its excellent transparency, rigidity, chemical and high temperature resistance make this material the most similar one to glass.

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**POLYETHYLENE OR POLYTHENE (PE)**  
*POLIETILENE O POLITENE (PE)*

It has good insulating properties and it is light, practically unbreakable, with an optimal dimensional and chemical stability. Solvent resistant, at ambient temperature no substance is able to melt PE. PE is not autoclavable.

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**POLYPROPYLENE (PP)**  
*POLIPROPILENE O MOPLEN (PP)*

It's an extremely versatile polymer, light, translucent, unbreakable, with an optimal mechanical resistance, but at the same time it's flexible. It has a good heat resistance (it can be repeatedly autoclaved at +121 °C) and an optimal dimensional stability. Optimal chemical resistance to solvents (at ambient temperature no solvent is able to melt it), to moisture, oils.

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**POLYCARBONATE (PC)**  
*POLICARBONATO (PC)*

It is as transparent as glass and very bright. In the thermal gap between -150 °C and +135 °C it resists to impact and it has high rigidity and solidity. It has a high dimensional stability as well as a remarkable resistance to rupture under stress. It has optimal optical properties and high UV resistance. PC resists to weak acids, mineral acids, halogens, aliphatic hydrocarbons, gasoline, fats, oils, water under +70 °C and to alcohols except for methyl alcohol. It can be sterilized in autoclave at + 121 °C for 20 minutes.

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It's the most resistant and rigid thermoplastic material and it has an optimal dimensional stability since it absorbs little moisture. It has a good mechanical resistance and an elevated hardness that gives an excellent abrasion and wear resistance. It has an excellent chemical resistance to organic solvents, weak acids, alkali, carburants, mineral oils, decalcifying solutions and solvents.

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It's commonly known as Teflon®. It's smooth, a little rigid and solid, tough and flexible at low temperatures. Its advantages are withstanding wide temperature ranges (it gets weak at temperature lower than -260 °C), universal chemical resistance, insoluble in all solvents below + 300 °C, resistant to atmospheric agents, exceptional heat stability (nonflammable) and anti-adhesive properties.

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**SILICONE RUBBER (SI)**  
**GOMMA SILICONICA (SI)**

Normally silicone rubber is very flexible, resists to ageing, chemical attacks and oxidation, is an excellent electrical insulator and has optimal non-stick properties. It resists to very high temperatures, up to +180 °C for long lasting use and up to +250 °C for short lasting use in warm and dry ambient. SI products are not autoclavable because steam attacks rubber starting from a temperature of +100 °C. SI is biocompatible, inert, neither toxic nor inflammable.

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All our plastic Devices for laboratory are easy to clean, as they present washable surfaces, resistant both to high temperatures and to the most commonly used detergents. In case of light contaminations we suggest the use of neutral detergents (pH 7.0), while for heavy contaminations we suggest the use of alkaline detergents (pH up to 12.0). Never use any alkaline detergents on Devices manufactured in polycarbonate (PC) or polystyrene (PS). In order to accelerate decontamination process, rise detergent's concentration and solution's temperature. Agitation and rubbing speed up decontamination. In any case, you are requested to check with your detergents supplier their compatibility with plastic materials of the Devices. Always avoid the use of powders or abrasive sponges when washing and cleaning any Devices. In case the Devices get contaminated by dangerous, infected or toxic products, we suggest you to proceed first with a decontamination and washing process and an appropriate sterilization. Autoclaving is the favourite sterilization method, but other chemical and dry heat methods can be employed as long as they are compatible with the Device's plastic. Waste disposal: used disposables have to be considered as hazardous waste and have to be disposed according to laws in force.

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Our products, according to European Norms in force and to destination of use, can be classified as follows:

- **In Vitro Medical-Diagnostics Devices (IVD):** reference is made to Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998. By IVD it's meant all Devices used for in vitro diagnostic tests of samples coming from human beings (veterinary samples excluded), included blood and donated tissues (it's not allowed using IVD directly on the patient). Notified Body is not necessary and Declaration of Conformity will be issued directly by the Manufacturer (except for Devices whose purposes are self diagnostic tests where the Notified Body has to be applied).
- **General Laboratory Use:** according to Directive 98/79/EC, those products that are not destined from the Manufacturer to contain directly a human sample for diagnostic tests. They are not CE marked.
- **Medical Devices (MD):** reference is made to Directive 93/42/EC of the European Council of 14 June 1993, Directive 2007/47/CE of the European Parliament and of the Council of 5 September 2007 and subsequent amendments and integrations. These Devices are destined to be used directly on humans. MD are divided into classes as follows: Class I) Non invasive. Class IIa and IIb) Invasive for natural orifices (IIa) and invasive of surgical type (IIb). Class III) Long time invasive or implantable.
- **Individual Protection Devices (IPD):** reference is made to Directive 89/656/EC (use destination) and 89/686/EC (planning). Moreover IPD have to conform to specific harmonized Norms. By IPD are meant all equipments that a worker has to wear in order to protect himself from any risks that could damage his safety or health during working hours. IPD can be divided into three categories, according to the risk type.
- **Medical-Surgical Instruments (MSI) (Only for Italian Law):** subject to D.P.R. n. 392/98, they are those products used in medicine, zootechnics and agriculture (for example disinfectants, bactericides or germicides) that can be produced and introduced to the market after authorization of Ministry of Health that grants their efficacy, reliability and safety during use.

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#### MICROBIOLOGICAL STATUS STATO MICROBIOLOGICO

Our products can be sold with the following microbiological status:

- **Non sterile:** Devices that are neither produced in contamination controlled conditions nor subjected to a following sterilization process, therefore they cannot be defined as "sterile" or "aseptic". Anyway Bioburden (medium microbial charge) of Devices is constantly controlled in order to verify contamination level of Devices to guarantee the Device quality.
- **Aseptic:** Devices produced and packed under controlled contamination conditions. For these Devices it's granted a SAL (Sterility Assurance Level) of 10<sup>-3</sup> for integral packaging. Every lot of material is subject to specific tests in order to verify microbiological status and afterwards deliberated.
- **Sterile:** by "sterile" it's meant a Device that is devoid of vital organism. Theoretically only a sterilization process that lasts forever could assure that all microorganisms get destroyed. A Device is defined as "sterile" when it's possible to assure or certify, on statistic basis, that the probability that it results non sterile is equal or lower than a chance of 1 in 1.000.000 i.e. SAL (Sterility Assurance Level) of 10<sup>-6</sup>, that is the one we guarantee for our "sterile" products. All our sterilization processes, carried out by irradiations of ionizing radiations (Gamma and Beta rays) or by Ethylene Oxide, are validated and constantly monitored in accordance with ISO 11137, ISO 11737, EN 552, F.U.I:XI and F.EU:IV Norms. The packaging of sterile Devices is validated and constantly monitored in order to guarantee the microbiological status for the indicated validity period of the integral packaging.
- **Apyrogenic or non-Pyrogenic:** Devices where absence or presence within detection limits of endotoxins by Gram-negative bacteria is granted. For Devices declared as Apyrogenic a limit of ≤ 20 EU/piece (EU=Endotoxin Unit) FDA Endotoxin on integral packaging is granted. Endotoxins presence is determined by LAL tests.

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**LEGAL DISCLAIMER - End Users** are responsible for the Devices they select and shall verify both suitability and completion of technical information provided, as products must comply in material, shape and dimension to their intended use. Nuova Aptaca S.r.l. has faculty to modify, without any previous notice, technical and aesthetic features of the articles described in this present catalogue.

All data and technical information provided are based on knowledges available at the time of the catalogue printing and to be intended for indicative purpose only and do not imply any responsibilities from Nuova Aptaca S.r.l. side. All information present on this catalogue exclusively belong to Nuova Aptaca S.r.l. and is strictly prohibited any partial or total reproduction. Aptaca®, CliniSwab®, VacuAptaca® are registered trade mark of Nuova Aptaca S.r.l., all right reserved.

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© Copyright Aptaca 2018. Tutti i diritti riservati.