
TECHNICAL MANUAL
MANUALE TECNICO

.25

APETACCA

QUALITY ASSURANCE AND COMPLIANCE OF APTACA DEVICES

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All Aptaca devices are designed, manufactured, verified, and approved according to rigorous quality procedures, with controls carried out using specific methodologies. Aptaca S.p.A. operates under a Quality System certified in compliance with **UNI EN ISO 9001, UNI CEI EN ISO 13485**, and other relevant certifications. Our devices meet the essential safety and compliance requirements set forth by applicable standards, directives, and regulations.

Our quality control processes include functionality tests, dimensional tests, physical tests, mechanical strength tests, leakage tests, microbiological tests, sterility tests, and regular inspections to ensure maximum safety and reliability of our devices. Each device is assigned to a **production batch** that is fully traceable, from raw material to finished product.

We offer technical support and customer service through various channels, ensuring a fast response times. For each device placed on the market, a **Certificate of Analysis, Quality, and Compliance (CoA)** is available, which can be requested from our Quality and Regulatory Office (qra@aptaca.com - regulatory@aptaca.com) or downloaded directly from our website www.aptaca.com upon registration. A wide range of **technical documents, guides, and educational resources** are also available to support the safe and effective use of our devices.

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MICROBIOLOGICAL AND MOLECULAR PURITY STATUS

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Our devices can be supplied with the following microbiological and molecular purity levels:

- **Non-sterile:** Devices neither produced under controlled contamination conditions nor subjected to subsequent sterilization processes, and therefore cannot be defined as "sterile" or "aseptic." The bioburden (average microbial load) of the devices is constantly monitored through periodic dose audits according to the standards of the current European Pharmacopoeia edition and ISO 11737-1, to ensure the quality of the device.
 - **Aseptic:** Devices produced and packaged under controlled contamination conditions. For these devices, a SAL (Sterility Assurance Level) of 10^{-3} is guaranteed with intact packaging. Each batch of material is subjected to specific tests to verify its microbiological status and subsequently released.
 - **Sterile:** A "sterile" device is free of viable organisms. A SAL of 10^{-6} means that the probability that a device is not sterile is less than 1 in 1,000,000. All our sterilization processes, carried out by irradiation with ionizing radiation (Beta or Gamma rays) or by Ethylene Oxide, are validated and constantly monitored. Beta ray sterilization processes are validated according to UNI EN ISO 11137, UNI EN ISO 11737, UNI EN 556, UNI EN 868, UNI EN ISO 11607, and periodic revalidation is performed according to ISO 11137-1, ISO 11137-2, ISO/TS 13004. Ethylene Oxide sterilization processes are validated, and revalidated periodically, according to ISO 11135, ISO 10993-7, and ISO 14937. The packaging of sterile devices is also validated and monitored in accordance with ISO 11607 to ensure microbiological integrity throughout the validity period with intact packaging.
 - **Pyrogen free or Non-pyrogenicity:** Devices guaranteed to be free or within detectable limits of endotoxins from Gram-negative bacteria. For devices declared apyrogenic, a limit of ≤ 20 EU/piece (Endotoxin Units) with intact packaging is guaranteed. The detection of endotoxins is carried out using the LAL Test (Limulus Amebocyte Lysate Test), a method that uses amebocytes from *Limulus polyphemus* (horseshoe crab) to detect the presence of bacterial endotoxins.
- Note:** Sterility and apyrogenicity are distinct concepts. Sterility refers to the absence of viable organisms, while apyrogenicity indicates the absence of endotoxins, which are substances capable of causing fever. A sterile product is not necessarily apyrogenic, as sterilization eliminates microorganisms but does not always remove endotoxins. Similarly, an apyrogenic product is not necessarily sterile.
- **Human DNA free:** Devices guaranteed to be free of human DNA contamination. Verification is performed using the Real Time PCR method, with a confidence level of approximately 95%. These devices are ideal for applications where the presence of contaminating human DNA could interfere with results, such as in molecular biology analyses.
 - **DNase and RNase free:** Devices guaranteed to be free of DNase and RNase contamination, enzymes that degrade DNA and RNA respectively. Verification is performed using fluorimetric assay, with a confidence level of approximately 95%. These devices are essential for molecular biology experiments and other applications where maintaining the integrity of DNA and RNA is crucial.
 - **ATP free:** Devices guaranteed to be free of ATP (adenosine triphosphate) contamination. Verification is performed using bioluminescence tests, sensitive to the presence of very low levels of ATP. These devices are used in applications where the presence of ATP could interfere with biochemical analysis results.
 - **PCR inhibitors free:** Devices guaranteed to be free of PCR inhibitors. Verification is performed using PCR amplification tests to confirm the absence of inhibitors that could interfere with the efficiency of the reaction. These devices are crucial for ensuring the effectiveness of PCR (Polymerase Chain Reaction) reactions, avoiding false negative results.

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Aptaca devices are manufactured with top-quality, medical-grade raw materials, ensuring quality and safety standards in line with the strictest national, European, and international standards. Production follows GMP (Good Manufacturing Practices).

In general, our devices are:

- **Suitable for food contact:** Made with materials that do not contain substances subject to restrictions according to Regulation 1895/2005/EC, Regulation 10/2011 (migration limits), Regulation (EU) 2020/1245, Regulation (EU) 2019/1338, and subsequent amendments. Compliance with overall and specific migration limits (where applicable). The raw materials used comply with US FDA Food Contact Chapter 21 CFR 177. This statement is supported by analytical tests conducted at qualified laboratories and information provided by the raw material suppliers.
- **Biologically inert and non-cytotoxic:** The devices are made from materials that do not release toxic or reactive substances that could harm human cells or tissues. Biocompatibility is verified according to ISO 10993-1 and ISO 10993-5 standards, which cover biological evaluation of medical devices and in vitro cytotoxicity tests. Tests are conducted in accredited laboratories to ensure that materials do not cause adverse reactions when in contact with the human body.
- **Hemocompatible:** Non-hemolytic, tested by accredited external laboratories according to the requirements of ISO 10993-4 and ASTM F756-17 standards.
- **REACH compliant:** In accordance with Regulation (EC) No. 1907/2006, they do not contain SVHC substances listed in the latest "Candidate list" (<https://echa.europa.eu/en/substances-restricted-under-reach>).
- **California Proposition 65 compliant:** They do not contain substances listed in the latest edition of "The Proposition 65 list".
- **Latex-free:** Manufactured without natural rubber, latex, or synthetic rubbers containing natural rubber (excluding latex items). Verification is performed using FT-IR spectra in ATR (attenuated total reflectance) mode to detect the presence of natural rubber latex with a sensitivity of 0.5%. This statement is based on tests conducted by qualified external laboratories and information provided by the raw material suppliers.
- **Free from heavy metals:** Compliant with Directive 2011/65/EU (RoHS II) and the amendments of Directive 2015/863/EU.
- **Phthalate and sulfate-free:** Compliant with Delegated Directive (EU) 2015/863 and Directive 2009/48/EC.
- **BSE/TSE-free:** Do not use animal-origin additives as defined by Regulation (EU) 10/2011 and are free from Specific Risk Materials (SRM) related to Bovine Spongiform Encephalopathy (BSE) and Transmissible Spongiform Encephalopathy (TSE).
- **Bisphenol-free:** The raw materials used are declared free from Bisphenol A, B, F, and S.
- **Halal/Kosher:** Ethanol (CAS No. 64-17-5) and animal-origin additives are not used in the production.
- **Compliant with cosmetic regulations:** Do not contain prohibited ingredients as indicated in Annexes II and III of Regulation (EC) 1223/2009 on cosmetic products.

For more information on legislation, applicable regulations, and specific declarations of conformity for each device, please contact our Quality and Regulatory Office (qra@aptaca.com - regulatory@aptaca.com).

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Each sales unit of our devices is clearly and correctly identified through a detailed and easy-to-read labeling system. We use **symbols in compliance with the UNI CEI EN ISO 15223-1 standard** (Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied), ensuring complete traceability, including through GS1 EAN128 standard barcodes, and safe use of the device. The information provided with the devices complies with EN 1041 (Information to be supplied by the manufacturer of medical devices) and ISO 20417 (Medical devices - Information to be supplied by the manufacturer) standards. Where applicable, the label includes the **UDI code** and QR code to access the **eIFU (electronic instructions for use)**, in accordance with Regulation (EU) 2021/2226.

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DECONTAMINATION, WASHING AND WASTE OF DISPOSABLE LABORATORY DEVICES

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All plastic laboratory devices are easy to wash, as they have wettable surfaces, resistant to both high temperatures and common chemical detergents used. For light contamination, washing with neutral detergents (pH 7.0) is recommended, while for more substantial contamination, the use of alkaline detergents (pH up to 12.0) is advised. Never use alkaline detergents with polycarbonate (PC) or polystyrene (PS) products. To accelerate decontamination, increase the concentration of the detergent and the temperature of the solution. Agitation and scrubbing can also speed up decontamination. In any case, check with the detergent supplier for compatibility with the plastic materials of the device. For washing and cleaning devices, always avoid using powders or abrasive sponges. In case of contamination of devices with hazardous, infectious, or toxic products, it is advisable to proceed with decontamination, washing, and adequate sterilization. Autoclaving is the preferred method of sterilization, but chemical or dry heat methods can also be used, provided they are compatible with the plastic used. Avoid cross-contamination by using dedicated equipment.

Disposal: Used devices are considered potentially infectious waste and must be disposed of according to current regulations. Unused devices, on the other hand, can be disposed of through regular recycling channels, following local guidelines for plastic materials. Ensure that single-use devices are disposed of correctly and not reused.

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GENERAL OBLIGATIONS OF DISTRIBUTORS

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In accordance with **Regulation (EU) 2017/745 (MDR)** and **Regulation (EU) 2017/746 (IVDR)** of the European Parliament and of the Council, a Distributor is any natural or legal person in the supply chain, other than the manufacturer or the importer, who makes a device available on the market until the point of putting it into service. The activities of distributors include procuring, holding, and supplying devices. The Distributor must comply with all the obligations set out in Article 14 of the aforementioned Regulations.

Distributors ensure that, for the period during which a device is under their responsibility, the **storage and/or transport conditions** comply with the conditions set by the manufacturer. Additionally, it is the responsibility of the Distributor to maintain the **traceability of devices** in accordance with the aforementioned Regulations and ISO 13485 standards. European Commission guidelines, such as MDCG 2021-19, provide further details on traceability requirements, including the integration of UDI into the quality management system.

Furthermore, Distributors must:

- Ensure that their personnel are adequately trained and competent to handle medical and in vitro diagnostic devices in accordance with applicable regulations.
- Implement and maintain a quality management system, preferably in compliance with ISO 9001 or ISO 13485 standards.
- Have procedures in place for device recalls, as required by Article 14 of the MDR/IVDR.
- Continuously monitor the market and report any incidents or issues found with the devices, in accordance with Article 25 of the MDR/IVDR.
- Be prepared for audits and inspections by competent authorities and notified bodies, as required by Article 93 of the MDR and Article 88 of the IVDR.

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Our devices, according to current EU regulations and their intended use, can be classified as follows:

IN VITRO DIAGNOSTIC MEDICAL DEVICES (IVD): The reference regulation is Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017, which has repealed Directive 98/79/EC and Commission Decision 2010/227/UE. This regulation establishes the rules for placing in vitro diagnostic medical devices on the market, making them available on the market, or putting them into service for human use and accessories for such devices in the Union. IVD refers to any device, used alone or in combination, intended by the manufacturer to be used in vitro for the examination of specimens derived from the human body, including blood and tissue donations, solely or principally for the purpose of providing information: 1) concerning a physiological or pathological process or state; 2) concerning a congenital physical or mental impairment; 3) concerning the predisposition to a medical condition or a disease; 4) to determine the safety and compatibility with potential recipients; 5) to predict treatment responses or reactions; 6) to define or monitor therapeutic measures.

IVDs are divided into classes A, B, C, and D, based on the intended use of the devices and the risks involved. Invasive products, intended for sample collection or products that come into direct contact with the human body to obtain a sample, are not considered IVDs.

Aptaca's IVDs are classified as Class A, further divided into:

- **Class A non-sterile:** These devices do not require the intervention of a Notified Body for conformity assessment. The EU declaration of conformity is based on the manufacturer's self-certification. Manufacturers must ensure that the devices meet general safety and performance requirements, prepare technical documentation, and declare conformity by drafting and maintaining the EU declaration of conformity.
- **Class A sterile:** These devices require the intervention of a Notified Body for conformity assessment concerning sterility aspects. Conformity is certified through a CE certificate issued by the Notified Body. The manufacturer must follow the conformity assessment procedure based on the quality management system and control of the sterilization process (Annex IX of Regulation (EU) 2017/746).

Classification Rule V (Annex VIII, Rule 5). Class A devices are those that present low risk and include: 1) Reagents and reagent products, including related calibrators and control materials; 2) Instruments, software, materials, and other articles intended by the manufacturer to be used specifically to provide information resulting from in vitro examination of specimens derived from the human body; 3) Sample containers.

DEVICES FOR GENERAL LABORATORY USE: Regulations (EU) 2017/745 (Article 1, paragraph 5(b)) and (EU) 2017/746 (Article 1, paragraph 3) do not apply to products intended for general laboratory use or products exclusively for research, unless, by their characteristics, they are specifically intended by the manufacturer for in vitro diagnostic examinations. These products do not have CE marking.

MEDICAL DEVICES (MD): The reference regulation is Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017, which amended Directive 2001/83/EC, Regulation (EC) No 178/2002, and Regulation (EC) No 1223/2009, and repealed Council Directives 90/385/EEC and 93/42/EEC. This regulation establishes the rules for placing medical devices for human use and their accessories on the market, making them available on the market, or putting them into service in the Union.

Medical Device refers to any instrument, apparatus, appliance, software, implant, reagent, material, or other article, intended by the manufacturer to be used for human beings, alone or in combination, for one or more of the following specific medical purposes: 1) diagnosis, prevention, monitoring, prediction, prognosis, treatment, or alleviation of diseases; 2) diagnosis, monitoring, treatment, alleviation, or compensation for an injury or disability; 3) investigation, replacement, or modification of the anatomy or of a physiological or pathological process or state; 4) providing information through the in vitro examination of specimens derived from the human body, including blood and tissue donations.

A medical device does not achieve its principal intended action in or on the human body by pharmacological, immunological, or metabolic means, but its function may be assisted by such means. Medical devices are classified into classes I, IIa, IIb, and III, based on the intended use of the devices and the risks they entail.

Aptaca's medical devices are classified into the following categories:

- **Class I non-sterile:** These devices do not require the intervention of a Notified Body for conformity assessment. The EU declaration of conformity is based on the manufacturer's self-certification. Manufacturers must ensure that the devices meet general safety and performance requirements, prepare technical documentation, and declare conformity by drafting and maintaining the EU declaration of conformity.
- **Class I sterile:** These devices require the intervention of a Notified Body for conformity assessment concerning sterility aspects. Conformity is certified through a CE certificate issued by the Notified Body. The manufacturer must follow the conformity assessment procedure based on the quality management system and control of the sterilization process (Annex IX of Regulation (EU) 2017/745).
- **Class IIa sterile:** These devices require the intervention of a Notified Body for conformity assessment. Conformity is certified through a CE certificate issued by the Notified Body. The manufacturer must follow the conformity assessment procedure based on the quality management system and control of the sterilization and safety process (Annex IX of Regulation (EU) 2017/745).

PERSONAL PROTECTIVE EQUIPMENT (PPE): The reference regulation is Regulation (EU) 2016/425 of the European Parliament and of the Council of 9 March 2016 on personal protective equipment, repealing Council Directive 89/686/EEC. PPE is designed to protect against health and safety risks and is divided into three categories based on the type of risk. They must comply with specific harmonized standards.

BIOCIDES (BPR): The reference regulation is Regulation (EU) No 528/2012 of the European Parliament and of the Council of 22 May 2012 concerning the making available on the market and use of biocidal products. A biocide refers to any substance or mixture intended to destroy or make harmless dangerous organisms. Production and marketing require authorization from the Ministry of Health, which ensures their effectiveness and safety.

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GENERAL FEATURES OF PLASTIC MATERIALS

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Aptaca devices are made from virgin, high-quality plastics specifically selected for medical use. We ensure that all materials used meet the strictest quality and safety standards at the national, European, and international levels.

For more detailed information about the technical specifications of the plastics used, please contact our Quality and Regulatory Office (qra@aptaca.com - regulatory@aptaca.com).

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ACRYLONITRILE BUTADIENE STYRENE (ABS)

ACRILONITRILE BUTADIENE STIRENE (ABS)

Compared to crystal polystyrene, that belongs to the same family, it's a more technical material and has a chemical resistance and temperature stability undoubtedly better. Products made of ABS are light and have a good rigidity grade. The most important mechanical properties are the high impact resistance and toughness also at low temperatures, hardness and scratch resistance. ABS is not autoclavable.

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POLYSTYRENE (PS)

POLISTIROLO O POLISTIRENE (PS)

It's a particularly transparent and bright material. It's rigid and hard but also rather fragile and subject to breakages. It is insensitive to moisture, steady at saline solutions, alkali, non oxidizing acids, aqueous chemical agents, but it has a very limited resistance to solvents. It's necessary to take into consideration that PS products, at room temperature, are particularly rigid and therefore they can crack or break in case they fall down from a laboratory bench or from higher heights. PS is not autoclavable.

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POLYMETHYL METHACRYLATE (PMMA)

POLIMETILMETACRILATO (PMMA)

PMMA is hard, rigid, quite fragile but extremely abrasion resistant (scratch resistant), it has high flexibility, it is stable to UV rays, but easily inflammable. PMMA has a good impact strength, higher than glass, but significantly lower than polycarbonate. Thanks to its optical properties, this plastic material is the most similar to glass. PMMA is not autoclavable.

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POLYMETHYLPENTENE (PMP-TPX®)
POLIMETILPENTENE (PMP - TPX®)

It's a very light resin, transparent (90% of transmittance of visual light) and very hard, with good mechanical properties. PMP resists to concentrated sulfuric acid, acetone, ethanol, mineral acids, strong and weak acids, strong and weak alkaline solutions, inorganic salts, aldehydes, alcohols, detergents, oils, fats and boiling water. It's steady to ketones with a limited extent. It has an optimal resistance to high temperatures (it can be repeatedly autoclaved both at +121 °C and at +150 °C). Its excellent transparency, rigidity, chemical and high temperature resistance make this material the most similar to glass.

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POLYETHYLENE OR POLYTHENE (PE)
POLIETILENE O POLITENE (PE)

It has good insulating properties and it is light, practically unbreakable, with an optimal dimensional and chemical stability. Solvent resistant, at ambient temperature no substance is able to melt PE. PE is not autoclavable.

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POLYPROPYLENE (PP)
POLIPROPILENE (PP)

It's an extremely versatile polymer, light, translucent, unbreakable, with an optimal mechanical resistance, but at the same time it's flexible. It has a good heat resistance (it can be repeatedly autoclaved at +121 °C) and an optimal dimensional stability. Optimal chemical resistance to solvents (at ambient temperature no solvent is able to melt it), to moisture, oils.

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POLYCARBONATE (PC)

POLICARBONATO (PC)

It is as transparent as glass and very bright. In the thermal gap between -135 °C and +121 °C it resists to impact and it has high rigidity and solidity. It has a high dimensional stability as well as a remarkable resistance to rupture under stress. It has optimal optical properties and high UV resistance. PC resists to weak acids, mineral acids, halogens, aliphatic hydrocarbons, gasoline, fats, oils, water under +70 °C and to alcohols except for methilic alcohol. It can be sterilized in autoclave at +121 °C for 20 minutes.

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POLYOXYMETHYLENE OR ACETAL RESIN (POM)

POLIOSSIMETILENE O RESINA ACETALICA (POM)

It's the most resistant and rigid thermoplastic material and it has an optimal dimensional stability since it absorbs little moisture. It has a good mechanical resistance and an elevated hardness that gives an excellent abrasion and wear resistance. It has an excellent chemical resistance to organic solvents, weak acids, alkali, carburants, mineral oils, decalcifying solutions and solvents.

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POLYTETRAFLUOROETHYLENE OR TEFLON® (PTFE)

POLITETRAFLUOROETILENE O TEFLON® (PTFE)

It's commonly known as Teflon®. It's smooth, a little rigid and solid, resistant and flexible at low temperatures. Its advantages are withstanding wide temperature ranges (it gets weak at temperature lower than -260 °C), universal chemical resistance, insoluble in all solvents below +300 °C, resistant to atmospheric agents, exceptional heat stability (not inflammable) and anti-adhesive properties.

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SILICONE RUBBER (SI)

GOMMA SILICONICA (SI)

Normally silicone rubber is very flexible, resists to ageing, chemical attacks and oxidation, is an excellent electrical insulator and has optimal non-stick properties. It resists to very high temperatures, up to +180 °C for long lasting use and up to +250 °C for short lasting use in warm and dry ambient. SI products are not autoclavable because steam attacks rubber starting from a temperature of +100 °C. SI is biocompatible, inert, neither toxic nor inflammable.

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POLYETHYLENE TEREPHTHALATE (PET)
POLIETILENE TEREFALATO (PET)

PET is a widely used plastic material known for its excellent transparency, strength, and dimensional stability. It has a smooth surface with great optical properties, making it ideal for applications where good content visibility is required. It has high mechanical strength, even under intensive use. PET is chemically resistant to many acids, bases, and saline solutions, but it is sensitive to strong solvents like alcohols and ketones. PET provides good rigidity and durability, making it suitable for producing vacuum tubes, ensuring sterility and vacuum retention, as well as bottles. PET is not autoclavable and cannot withstand temperatures above 70°C without deformation.

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LEGAL DISCLAIMER
NOTE LEGALI

The end user is responsible for choosing devices that are suitable in terms of material, shape, and size, according to their intended use, and for ensuring the conformity and completeness of technical information. Aptaca S.p.A. reserves the right to modify, without prior notice, the technical and aesthetic characteristics of the items listed in the catalogue. All data and technical information provided are based on the knowledge available to us at the date of printing and are intended for indicative purposes only, implying no responsibility on the part of Aptaca S.p.A. The end user is also responsible for ensuring that the use of the devices complies with all applicable laws, regulations, and standards in their country. Aptaca S.p.A. will not be held liable for any improper use of the devices, including but not limited to, uses that do not conform to the provided instructions or uses in unsuitable environments.

All information in this catalogue is the exclusive property of Aptaca S.p.A. and reproduction, in whole or in part, without prior written consent is strictly prohibited. Aptaca®, CliniSwab®, VacuAptaca® are registered trademarks of Aptaca S.p.A., all rights reserved.
© Copyright Aptaca 2024, all rights reserved.

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+= Completely resistant / Resistence ~=Resistant with expection / Mediamente resistente - = Not resistant / Non resistente

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© Copyright Aptaca 2024, tutti i diritti riservati.