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### Bio Transport Kit UN3373

*Complete system with a rigid secondary packaging and a rigid outer packaging suitable for the triple packing transport system of 4 test tubes containing Biological Substance Category B, UN 3373. An example of "Biological Substance Category B" is swabs sent for the Covid-19 (Coronavirus) test.*

*System compliant with packaging instructions P650 ADR (for shipments by land), P650 RID (for shipments by rail), P650 ADN (for shipments by inland waterway), P650 IMDG (for shipments by sea) and P.I. 650 IATA (for shipments by air) and as prescribed by the World Health Organization "Guidance on regulations for the transport of infectious substances 2019 - 2020". Compliant with the requests of the Ministry of Health regarding the containers to be used for the transport and shipment of swabs for Coronavirus tests (Covid-19). Possibility to insert cooling element cod. 6801 (not included in the kit).*

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Infectious substances are divided into two categories:

**Category A:** An infectious substance transported in such a state that, when exposed to the same, can produce a permanent disability, it is potentially fatal or can cause a fatal illness in humans or animals. Code UN 2814 will be assigned to infectious substances that cause diseases among humans or both humans and animals. UN 2900 code will be assigned to infectious substances which cause disease only in animals.

**Category B:** any infectious substance which does not fulfill the criteria of inclusion in Category A including diagnostic specimens and clinical. Code UN 3373 will be assigned to infectious substances of Category B to the. An example of "Category B Biological Substance" is swabs sent for the Covid-19 (Coronavirus) test.

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### PACKAGINGS FOR BIOLOGICAL SUBSTANCE, CATEGORY B (UN3373)

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*For the shipment of samples such as blood or tissues swabs (example for Covid-19), when they fall under the UN classification "Biological Substance, Category B", they need an appropriate packaging in accordance with the regulations for dangerous goods.*

*In fact, Category B Biological Substance requires packaging made up of a primary container, a secondary container and an external packaging.*

*The legislation on dangerous goods requires that these samples (blood, urine, tissues, etc...) are placed in a primary container that must be placed in a secondary packaging certified 95kPa from +40°C to +55°C as required by the ADR, RID, ADN, IMDG and ICAO / IATA regulations.*

*Everything must be placed in an external packaging that must be able to withstand the drop test from 1,2 meters high as required by the P650 ADR (for shipments by road), P650 RID (for shipments by rail), P650 ADN (for shipping by waterway), P650 IMDG (for shipments by sea) and P.I. 650 IATA (for shipments by air).*

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*Keep out of flame or heat sources which might damage the product*

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*Do not re-use: Disposable device*

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*Do not vary the intended purpose of the product*

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
*Store in dry place, Temperature range: min -10°C max +50°C*

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*Disposal: use appropriate personal protective equipment and act according to applicable regulations*

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- 1) Check that the primary container (test tube) is tightly closed. If in doubt about the seal, we recommend sealing with Parafilm® film;
- 2) Wrap the primary container in sheet of absorbent paper and insert into the secondary container. Close the secondary container by screwing the cap firmly;
- 3) Apply the SecurLabel tamper evident labels on the secondary containers identifying its content (optional);
- 4) Insert the absorbent foam support on the bottom of the tertiary container (outer packaging) in cardboard;
- 5) Insert the secondary containers in the appropriate housings of the absorbent sponge support;
- 6) Insert the accompanying documents inside the tertiary container (outer packaging);
- 7) Close the tertiary container (outer packaging) and seal it with the appropriate label contained in the kit;
- 8) Write the sender (Shipper) and recipient (Receiver) data on the external packaging.
- 9) The shipper (may also be known as the consigner or the sender):
  - ensures the correct classification, packaging, labelling and documentation of all infectious substances destined for transport;
  - ensures that packaging selected is suitable and compliant for the substances being shipped.

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*UN 3373, Biological Substance category B*

*Conformity to P650 ADR/RID*

- *internal pressure tests (hydraulic) in accordance with paragraph 6.1.5.5 of the ADR 2013 standard*
- *drop tests in accordance with paragraph 6.1.5.3 of the ADR 2013 standard*

**Product not subject to CE marking. For professional use only**

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Quantity (pcs): n. 5 kits suitable for the total transport of 20 primary containers (test tubes)

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MINIMUM SALEABLE QUANTITY

Misura esterna scatola (cm): 29,3 x 17,5 x 24  
External box dimensions (cm):

Peso (Kg): 0,95  
Weight (Kg):

Volume (m<sup>3</sup>): 0,012  
Volume (m<sup>3</sup>):

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**Aptaca S.p.A.** Regione Monforte, 30 - 14053 Canelli (Asti) Italy